# **Buntes Lern-Mikroskop** Nr./Art. 757943

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Datum \_\_\_\_\_

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## Warning notice for a safe and correct method of use

This microscope is not suitable for children under 8 years, since it contains small parts, which pose a choking hazard to young children. It is not a toy! It is pointed out to parents and teachers that the microscope is a scientific device. Cover slips, microscope slides, and sharp instruments belong to the microscope. Children are not allowed to work with these, unless supervised by adults, as there is a possible danger of getting hurt. Always follow the specified safety regulations! With the colourful learning microscope made by Betzold you have the option to directly view specimens on microscope slides as well as objects such as plants, coins or insects.

This is possible through the dual LED illumination, which enables you to light the object from below (transmitted light) as well as from above (reflected light). Batteries make you independent from electric supply lines, and enables you to use it outdoors. There are two different application options in one set of the colourful learning microscope, promising many hours of fun and learning.

## Maintenance, care, and storage of the microscope

**WARNING:** To avoid fire or short circuit, never expose the microscope to dampness or moisture.

- Turn the light off whenever you do not need it, in order to prolong the life of the batteries. Make sure, both lights are turned off, e.g. the switch is in the middle position.
- Clean the microscope exclusively with a soft cloth. Fingerprints and stains can be removed with a moistened cloth. When it is necessary to clean the lenses, try to do it with a cotton-wool ball, very slightly moistened with alcohol. Dry the lenses with a dry cotton-wool ball or with special cleaning paper.
- Store the microscope in a cool and dry place.
- Cover the microscope with the included dust cover, or put it back into its case whenever not in use.
- When you use it outside, make sure it is in an upright position, to prevent the eyepiece from falling out.
- Eyepiece and objective should be kept installed, to prevent dust getting into the tube.
- Store the unused microscope slides in the provided container. Always be careful with glass and sharp objects and keep them inaccessible for small children.
- When the light begins to weaken, replace the type AA (Mignon LR6) batteries. The LED lights have life time durability and should not have to be replaced.
- Normal batteries are not rechargeable! Replace used batteries with new ones.
- Do not mix old and new batteries.
- Do not keep used batteries in the battery holder for a longer period of time.
- Only use batteries of the same type.
- Please take care to match the right polarities.
- Do not expose batteries to fire or high temperatures!

## Notes for disposal

Please dispose of all packaging material immediately after unpacking in an environmentally-friendly fashion. Plastic wrappings represent a danger of suffocation for small children. Please dispose of the used microscope according to local ordinance regulations.

## Electrical waste is not to be disposed of in household waste.

If you cannot use the microscope any more, the consumer is legally bound to dispose of the used electrical devices separately from the household waste e.g. at a collection point of his municipality / city. This is indicated by the following symbol.

## Notes for disposal of batteries

Used batteries are not to be disposed of with household waste. Consumers are legally obligated to return batteries to the appropriate collecting points. They can also be taken to any place where they are sold. Used batteries contain valuable resources which can be recycled.

The sign with the crossed out waste container means: Batteries and rechargeable batteries are not to be disposed of in household waste.

The symbols below the waste container mean:

Pb: Battery contains lead

Cd: Battery contains cadmium

Hg: Battery contains mercury

## Guarantee

Over and above the legal warranty (and without restricting it) you receive 5 years complete guarantee. That means, proof of damaged or faulty goods at the time of purchase is not necessary. In case of guarantee claims, please contact place of purchase.

## Technical dates and delivery schedule

### a) Magnifications

The total magnification is calculated by multiplying (always by ten) the objective magnification with the size of the objective lens.

Objective	4x	10x	40x
Eyepiece (lower light)	40x	100x	400x
Stereo (light above)	40x	100x	not possible

Note: The 40x objective is not suitable for observations with reflected light. The eyepiece is then too close to the object; so there is not enough light for observation.

### b) Equipment

The colourful learning microscope is equipped as follows:

- 10x eyepiece magnification
- revolving nosepiece
- reflected LED light
- transmitted LED light
- focus knob for the LED illumination
- orifice plate
- optical glass lenses
- dimensions (H x W x D): 29,5 x 11,5 x 16,5 cm
- works with 3 type AA Mignon batteries (not included)

### c) Components of the colourful learning microscope

Carefully study the picture, and get yourself acquainted with the different components of your microscope.

### d) Accessories

- 5 empty slide holders
- pipette
- 4 prepared slide holders
- 2 small bottles of non poisonous colouring liquid
- test tubes
- 1 slide holder with dip
- petridish
- labels for slide holders
- tweezers
- cover glasses
- special cleaning paper for lenses
- skewer
- swab

**ATTENTION!** The set contains food colour and sharp instruments! Children are only allowed to do experiments and slide preparations when supervised by adults!

## Instructions

### a) Preparation

Please read this text carefully before using the microscope. Familiarise yourself with the individual components of the microscope.

#### Installing the batteries:

Detach the eyepiece from the lens tube and turn the microscope. Open the battery holder at the bottom and insert 3 batteries type AA (Mignon LR6) (not included) as indicated. Shut the battery holder and reinstall the eyepiece.

Always put your microscope on an even underground, such as a table.

The **colourful learning microscope** allows to illuminate the object from the top or from the bottom. Note that the switch at the back of the microscope serves for both, the transmitted as well as the reflected light, depending on the direction into which you turn it. The two lights cannot be used at the same time.

The lens tube and the top can be turned in such a way that they, if wanted, point to the other side of the microscope. If necessary, you can LIGHTLY turn the screw at the top base to make the turning of the top easier. For best results the top should only be moved when it is absolutely necessary.

### b) Work with transmitted light: Use of the bottom light source (for slide preparations on slide holders)

Put the stage into its lowest position. Always begin your observations with the 4x objective and increase the magnification gradually. Start with a prepared object holder and place it cautiously on the stage, the labelled side towards the top. Press the clamping device lightly in order to lift it up; to be able to fasten the object holder underneath. The part to be observed of the object should be placed on the opening in the stage. Use the switch for the bottom light so that it shines through the opening and through the object holder onto the stage. Look through the eyepiece and turn the focus knob SLOWLY, until the object can be seen clearly. Make sure the object holder does not touch the objective.

#### Hints for working with transmitted light:

- The orifice plate which sits directly under the stage has 6 different openings, which are used to control the light quantity. The turning of this disc enables you to change the intensity of the light. Try to experiment with the different possibilities of adjustment to get the best sight.
- After adjustment you can move the object holder to the left and to the right on the stage to observe different areas of the object. When you are adjusting the object holder or the disc, pay attention not to put any pressure, which could move the specimen onto the stage. In case this should happen all the same, simply readjust the focus knobs until you can see the object again clearly.
- After having looked at the 40x magnified specimen, (4x objective) you turn the objective turret to the 10x objective. (The object holder stays on the stage.) When SLOWLY turning the focus knob now you will be able to see your specimen magnified 100x. Always be very cautious and take care that the object holder does not touch the objective. Finally you can turn the objective turret to the 40x lens, which gives you a 400x magnification of your specimen. Whenever you increase, consider: The bigger the magnification, the nearer the lens must be to the object. The 40x objective almost seems to touch the lens; that is normal. Always move the focus knobs very slowly.
- When changing the object holder, put the stage into its lowest position. Lift the clips and remove the object holder. In order to fasten the new object holder, simply lift the clips again and push the object holder carefully underneath. Begin your observations again with the 4x objective.

**c) Work with reflected light: Use of the above light source  
(for solid objects and for the observation of surface structures)**

Switch the above light on to illuminate the stage. Just as with the work with the transmitted light, always begin with the 4x objective and increase the magnification accordingly.

**ATTENTION!** For the observation of solid objects only the 4x and the 10x objectives are suitable. The 40x magnification is unsuitable for the use with reflected light due to the necessary distance needed for the focus. For contrast reasons it might be necessary to place some specimen on a card or a small piece of paper before placing them on the stage. As you are using the light from above, it is not necessary that the opening in the stage stays open. Furthermore, you avoid your specimen falling through the hole.

This kind of observation makes it possible to magnify countless objects from daily life: e.g. small stones, leaves, insects, flowers, coins, stamps etc.

**Hints for working with reflected light:**

You will notice (for example with a little insect) that you can adjust the focus on different levels. This shows that it is not possible to focus the whole object with only one adjustment as it is possible with flat objects like stamps, where the whole object is on only one level.

## Suggestions for the observation of specimen

### CHILDREN ARE NOT ALLOWED TO DO EXPERIMENTS WITHOUT BEING SUPERVISED BY AN ADULT!

1. Place some grains of salt or sugar on the object holder.
2. Dissolve a teaspoon of salt/sugar in warm water in a test tube. Then use the pipette to place some drops on an object holder. Let them dry and observe the recrystallisation of the salt/sugar.
3. Shake some pollen of a blossom onto an empty object holder.
4. It is interesting to observe parts of feathers or any hair, that of domestic animals included. Try to compare the hair of different animals. You can fix the hair with a small piece of sticky tape onto the object holder.
5. Observe the fibres of different types of paper, or compare the print of a newspaper with that of a type writer or a laser printer.
6. Look at the colours in comics, or at strings and fibres of different materials.
7. Compare natural fibres like cotton wool or wool to synthetic ones like polyester.
8. Water out of a pond contains living organisms. Put a drop of it onto the object holder with the dip to observe them.
9. Very small insects can be looked at with little magnification.
10. Or try to look at parts of bigger beetles like legs or wings. Doing this, the needle is useful for the positioning. Always take insects with tweezers and wash your hands carefully after handling them.
11. Watch the mould on old bread or other food. Or create your own mould culture by putting a little piece of bread into a petri dish for some days. A little moistness accelerates the moulding process. (Here too, wash your hands carefully after work.)
12. Many objects seem to be transparent under a microscope. Usually, these objects are coloured to better recognize the cells. Iodine is a popular dye. You can also put your specimen in food colour for some minutes in order to dye them. Use the tweezers to handle the coloured objects. Let the coloured objects dry before putting them under the microscope. Keep in mind that the colour dyes all objects it gets in touch with, so watch out for your clothes and furniture.

## Suggestions for self-made specimen

Always start with clean object holders and cover glasses. In order to avoid fingerprints, always take the object holders exclusively by their sides. Place the cover glass cautiously using the tweezers.

**Special attention and supervision by adults is necessary** when the scalpel is used to cut objects. Whenever you do not use the scalpel, wrap it up and store it in a secure place.

### a) Production of short-term or “wet” specimens (e.g. tissue samples such as plant parts)

1. Use the scalpel to cut the part of the plant as thin as possible.
2. Place the specimen onto the object holder.
3. Use the pipette to place a drop of water onto the object holder.
4. Cautiously place the cover glass on top, so that the water spreads.
5. The cover glass should be placed as flat as possible to avoid bubbles.  
Caution: The thin cover glass will break when the pressure is too strong.
6. Carefully dab up eventually spilled water.
7. Now place the specimen carefully under the microscope (cover glass towards the top) and watch it then.
8. Use labels to mark your objects.

Maybe you prefer to exchange step two and three, in this case, you first put the drop of water onto the object holder and then the specimen, before adding the cover glass. Experiment with both methods to find the best one for your object. (Keep in mind that they are short-lived specimens).

#### Tip:

The object holders can be washed carefully, and then used again. The cover glasses are extremely thin and cannot be reused. Dispose of them carefully after use.

### b) Making simple long-term specimen

Proceed the same way as above, only try to use transparent glue instead of water to keep the object and the cover glass on the object holder. Let the glue dry before you look at the object. Use the labels to mark your specimen.

## Driving test for microscopy

Paint the following components of the microscope with the indicated colours. Label the components.

Lens barrel: grey


Stage: black

Specimen clamps: silver-grey

Illumination: yellow

Adjustment knob blue


Eyepiece: green

Solution on page 16.

# Driving test for microscopy

Label the accessories:

Solution: 1 Eyepiece, 2 Lens tube, 3 Reflected light, 4 Arm, 5 Orifice plate, 6 Focus knob, 7 Base, 8 Transmitted light, 9 Stage clip, 10 Stage, 11 Objective, 12 Objective turret, 13 Petri dish, 14 Colouring liquid, 15 Covering liquid, 16 Slide holders, 17 Test tubes, 18 Pipette, 19 Skewer, 20 Swab, 21 Tweezers

## Driving test for microscopy


Correctly assign the terms to the corresponding explanations:

Focus knob	1
Eyepiece	2
Plate	3
Base	4
Lens tube	5
Objective turret	6
Objektive	7
Stage	8
Object	9
Arm	10
Illumination	11
Stage clip	12

A	It is used to adjust the brightness of the picture.
B	The eyepiece is in this tube.
C	Here are the different objects. I can change the used objective by turning.
D	With this, the microscope stands safely.
E	I put the specimen which I would like to view on this.
F	I use this to generate a magnified picture of the specimen. On the socket it is stated to which size the picture is magnified.
G	It serves to carry the microscope.
H	It is used to illuminate the specimen.
I	It is used to focus on the picture.
J	They keep the specimen holders in place.
K	The lens through which I view into the microscope. On the socket it is stated how many times the picture is magnified.
L	The specimen which I would like to view.

Solution: 11, 2, K, 3, A, 4, D, 5, B, 6, C, 7, F, 8, E, 9, L, 10, G, 11, H, 12, J

## Driving licence for microscopy:

Name \_\_\_\_\_

can designate the parts of a microscope and operate a microscope according to the intended use.

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© Arnulf Betzold GmbH
Arnulf Betzold GmbH
Betzold Lernmedien GmbH