Schüler-Mikroskop PA 05

Anleitungsheft zu Nr./Art. 85905

Schüler-Mikroskop PA 05 LED

Anleitungsheft zu Nr./Art. 85910

Student Microscope PA 05

Instruction Manual for Nr./Art. 85905

Student Microscope PA 05 LED

Instruction Manual for Nr./Art. 85910

Microscope PA 05

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Microscope PA 05 LED

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Mikroskopu szkolnego PA 05

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Mikroskopu szkolnego PA 05 LED

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Betzold

Warnhinweise

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Warning notice for a safe and correct method of use

The microscope is designed for children ages 10 and up. It contains small parts, which pose a choking hazard to young children. It is not a toy! It is pointed out to parents and teachers that the microscope is a scientific device. Always follow the specified safety instructions. Only use the microscope with a power supply of 220 Volt. Connecting it to different power supplies damages the device.

Maintenance, care, and storage of the microscope

WARNING: To avoid fire or short circuit, never expose the microscope to dampness or moisture.

- Turn the light off whenever you do not need it, in order to prolong the life of lighting.
- Separate the microscope from the electric circuit before cleaning it.
- Clean the microscope exclusively with a soft cloth. Fingerprints and stains can be removed with a moistened cloth. When it is necessary to clean the lenses, try to do it with a cotton-wool ball, very slightly moistened with alcohol. Dry the lenses with a dry cotton-wool ball or with special cleaning paper.
- Store the microscope in a cool and dry place.
- Cover the microscope with the included dust cover, or put it back into its case whenever not in use.
- When you use it outside, make sure it is in an upright position, to prevent the eyepiece from falling out. Avoid dampness and moisture.
- Eyepiece and objective should be kept installed, to prevent dust getting into the tube.
- Always carry the microscope holding it at the arm.

Notes for disposal

Please dispose of all packaging material immediately after unpacking in an environmentally-friendly fashion. Plastic wrappings represent a danger of suffocation for small children. Please dispose of the used microscope according to local ordinance regulations.

Electrical waste is not to be disposed of in household waste.

The unusable microscope is not to be disposed of with household waste. Please dispose the used microscope according to local ordinance regulations. This is indicated by the following symbol.

Notes for disposal of batteries

Used batteries are not to be disposed of with household waste. Consumers are legally obligated to return batteries to the appropriate collecting points. They can also be taken to any place where they are sold. Used batteries contain valuable resources, which can be recycled. The sign with the crossed out waste container means, that batteries and rechargeable batteries are not to be disposed of in household waste.

The symbols below the waste container mean:

Pb: The battery contains lead

Cd: The battery contains cadmium

Hg: The battery contains mercury

Guarantee

Over and above the legal warranty (and without restricting it) you receive 5 years complete guarantee. That means, proof of damaged or faulty goods at the time of purchase is not necessary.

In case of guarantee claims, contact the place of purchase.

D: Arnulf Betzold GmbH; Ferdinand-Porsche-Str. 6; 73479 Ellwangen;
 Tel.: +49 7961 - 9000 - 0; Fax: +49 7961 - 9000 - 50; E-Mail: service@betzold.de; Internet: www.betzold.de

AT: Arnulf Betzold GmbH; Seebühel 1; 6233 Kramsach/Tirol;
 Tel.: +43 5337 - 644 - 50; Fax: +43 5337 - 644 - 59; E-Mail: service@betzold.at; Internet: www.betzold.at

CH: Betzold Lernmedien GmbH; Winkelriedstrasse 82; 8203 Schaffhausen;
 Tel.: +41 52 - 64480 - 90; Fax: +41 52 - 64480 - 95; E-Mail: service@betzold.ch; Internet: www.betzold.ch

Technical information and scope of delivery


a) Magnifications and specifications

Intensity of magnification = Objective x Eyepiece Lens

Achromatic objective	4x magnification	10x magnification	40x magnification with slide protector
Magnification with 10x wide field eyepiece	40x magnification	100x magnification	400x magnification

b) Components of Student-Microscope PA 05 and PA 05 LED

Study the picture below to become familiar with the different parts of your microscope

Equipment

- 1 Dust Cover
- 1 Spare Tungsten Bulb, 220/230 Volt, 20 Watt, for PA 05
- 1 Battery Recharger for PA 05 LED

Intended use of the microscope

The microscope PA 05 is particularly designed for teaching pupils microscopic viewing. It is a scientific device and not a toy.

Directions for use

a) Preparations: Setting up the microscope

Take the microscope out of its Styrofoam wrapping and put it onto a stable and level working surface.

b) Work Instructions

Illumination

The model PA 05 works with a tungsten bulb (220/230 Volt, 20 Watt) imbedded in the base. This bulb has an average operating time of 3000 working hours.

The model PA 05 LED is equipped with LED lights. The electricity for the lighting is stored in rechargeable batteries. The light of the LED bulbs hardly heats liquid specimens. This bulb has an average operating time of 100 000 working hours. The lighting is switched on at the back of the base. The rechargeable batteries are best recharged overnight.

Power supply

The model PA 05 LED can be connected to a normal power socket or used cable-free with the rechargeable batteries. One battery charge lasts for about 10 operating hours. The included recharger can only be connected to a power circuit of 220 Volt. The recharger is protected against overcharging. First connect the recharger to the microscope, and then connect it to the power circuit. Should the rechargeable batteries totally be discharged, please recharge them overnight. The batteries can only be recharged when the microscope is switched off. Should the rechargeable batteries be defective and need to be changed, please contact place of purchase.

The model PA 05 only operates with a power network of 220 Volt.

Placing the microscope slide

Place the specimen holder provided with a sample object onto the microscope stage (size 9,5 x 9,5 cm), and fasten it carefully with the help of the flexible spring clips.

Adjusting the image sharpness

Start by viewing your object with a four times magnification. Turn the knob for coarse adjusting until you receive a clear picture. Then use the knob for fine adjusting to achieve an optimal sharpness. Proceed the same way with every other object. Pay attention that the object does not touch the objective lens. In order to avoid damage to the objective lens or to the specimen preparation, the microscope is equipped with a specimen protection. This prevents the stage from being turned up too far and thus the objective lens pressing onto the object. The specimen protection is already adjusted correctly, but – if necessary – can be readjusted with the screw on the top side of the stage. Furthermore, the 40x objective is spring loaded, to avoid the specimen being pressed.

Regulating the light intensity

The luminance can be adjusted through the different openings of the aperture plate. Light influences the viewing quality, especially with high magnification.

Changing the objective

When you reached the optimal top view with the smallest magnification (4x), turn the objective lens piece to the next magnification (10x). The objective lens piece should click into position. Only a slight adjustment should be necessary, to receive a clear picture again. Should you increase the magnification further; the top view will be improved by adjusting the light.

After use

When having finished working with the microscope, turn off the current. Use the knob for coarse adjusting to move the object stage downwards. The device should be protected by a dust cover, or disassembled and put back into the original Styrofoam packing.

c) Exchanging the tungsten bulb of the microscope model PA 05

A tungsten bulb (220/230 Volt, 20 Watt) is required.

First pull out the power plug! Loosen the screws at the bottom of the microscope, and remove the base plate. Replace the used bulb by a new one. Pay attention – the bulb should not be touched with your fingers, as this might shorten the operating time. Replace the base plate.

Suggestions for the observation of specimen

Children should not experiment without adult supervision!

1. Place some grains of salt or sugar on the microscope slide.
2. Dissolve a teaspoon of salt/sugar in warm water in a test tube. Then use the pipette to place some drops on an object holder. Let them dry and observe the recrystallization of the salt/sugar.
3. Shake some pollen of a blossom onto an empty microscope slide.
4. It is interesting to observe parts of feathers, or any hair, that of domestic animals included. Try to compare the hair of different animals. You can fix the hair with a small piece of sticky tape onto the microscope slide.
5. Observe the fibres of different types of paper, or compare the print of a newspaper with that of a type-writer or a laser-printer.
6. Look at the colours in comics, or at strings and fibres of different materials.
7. Compare natural fibres like cotton wool or wool, to synthetic ones like polyester.
8. Water, out of a pond, contains living organisms. Put a drop of it onto the microscope slide with the dip to observe them.
9. Very small insects can be looked at with little magnification.
10. Or try to look at parts of bigger beetles, like legs or wings. Doing this, the needle is useful for the positioning. Always take insects with tweezers, and wash your hands carefully after handling them.
11. Watch the mould on old bread or other food. Or create your own mould culture, by putting a little piece of bread into a petri dish for some days. A little moistness accelerates the moulding process. (Here too, wash your hands carefully after work.)
12. Many objects seem to be transparent under the microscope. Usually these objects are dyed, to make the cells more visible. A current dye is iodine. You can also dip your tissue sample for some minutes into normal, diluted food colouring, in order to dye it. Use a pair of tweezers to transport the dyed objects. Bear in mind that the colour will dye any objects contacted. Therefore take care of clothing and furniture.

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