# Mikroskop M-TOP 600

Anleitungsheft zu Nr./Art. 86482

# Mikroskop M-TOP 600 LED

Anleitungsheft zu Nr./Art. 86483

## Microscope M-TOP 600

Instruction Manual for Nr./Art. 86482

## Microscope M-TOP 600 LED

Instruction Manual for Nr./Art. 86483

## Microscope M-TOP 600

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## Microscope M-TOP 600 LED

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## Mikroskopu M-TOP 600

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## Mikroskopu M-TOP 600 LED

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# Betzold

## Warnhinweise

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## Warning notice for a safe and correct method of use

This microscope is designed for children ages 9 and up. It contains small parts, which pose a choking hazard to young children. It is not a toy! It is pointed out to parents and teachers that the microscope is a scientific device. Before using this microscope, please read these instructions carefully in order to get to know the nature of the individual parts and to achieve the best possible results when working with the microscope. Keep these instructions. Always follow the specified safety instructions. Only use the microscope M-TOP 600 with a power supply of 220 Volt. Connecting it to different power supplies damages the device.

## Maintenance, care, and storage of the microscope

**WARNING:** To avoid fire or short circuit, never expose the microscope to dampness or moisture.

- Turn the light off whenever you do not need it, in order to prolong the life of lighting.
- Separate the microscope from the electric circuit before cleaning it.
- Clean the microscope exclusively with a soft cloth. Fingerprints and stains can be removed with a moistened cloth. When it is necessary to clean the lenses, try to do it with a cotton-wool ball, very slightly moistened with alcohol. Dry the lenses with a dry cotton-wool ball or with special cleaning paper.
- Store the microscope in a cool and dry place.
- Cover the microscope with the included dust cover, or put it back into its case whenever not in use.
- When you use it outside, make sure it is in an upright position, to prevent the eyepiece from falling out. Avoid dampness and moisture.
- Eyepiece and objective should be kept installed, to prevent dust getting into the tube.
- Always carry the microscope holding it at the arm.

## Notes for disposal

Please dispose of all packaging material immediately after unpacking in an environmentally-friendly fashion. Plastic wrappings represent a danger of suffocation for small children. Please dispose of the used microscope according to local ordinance regulations.

## Electrical waste is not to be disposed of in household waste.

The unusable microscope is not to be disposed of with household waste. Please dispose the used microscope according to local ordinance regulations. This is indicated by the following symbol.

## Notes for disposal of batteries

Used batteries are not to be disposed of with household waste. Consumers are legally obligated to return batteries to the appropriate collecting points. They can also be taken to any place where they are sold. Used batteries contain valuable resources, which can be recycled. The sign with the crossed out waste container means, that batteries and rechargeable batteries are not to be disposed of in household waste.

The symbols below the waste container mean:

Pb: The battery contains lead

Cd: The battery contains cadmium

Hg: The battery contains mercury

## Guarantee

Over and above the legal warranty (and without restricting it) you receive 5 years complete guarantee. That means, proof of damaged or faulty goods at the time of purchase is not necessary.

In case of guarantee claims, contact the following addresses:

D: Arnulf Betzold GmbH; Ferdinand-Porsche-Str. 6; 73479 Ellwangen;  
 Tel.: +49 7961 - 9000 - 0; Fax: +49 7961 - 9000 - 50; E-Mail: service@betzold.de; Internet: www.betzold.de

AT: Arnulf Betzold GmbH; Seebühel 1; 6233 Kramsach/Tirol;  
 Tel.: +43 5337 - 644 - 50; Fax: +43 5337 - 644 - 59; E-Mail: service@betzold.at; Internet: www.betzold.at

CH: Betzold Lernmedien GmbH; Winkelriedstrasse 82; 8203 Schaffhausen;  
 Tel.: +41 52 - 64480 - 90; Fax: +41 52 - 64480 - 95; E-Mail: service@betzold.ch; Internet: www.betzold.ch

## Technical Information and scope of delivery

### a) Magnifications and specifications

Objective: Achromatic objective

| Objective | Aperture (N.A.) |
|-----------|-----------------|
| 4-fach    | 0,10            |
| 10-fach   | 0,25            |
| 40-fach   | 0,65            |
| 60-fach   | 0,85            |

Intensity of magnification

|                        |      | Objective |       |       |
|------------------------|------|-----------|-------|-------|
| Lens & additional lens | 4 x  | 10 x      | 40 x  | 60 x  |
| 10 x working area      | 40 x | 100 x     | 400 x | 600 x |

Linsen

| Lens | Linear field of vision |
|------|------------------------|
| 10 x | 18 mm Diameter         |

Working area of stage


| Size         | X-direction moving distance | Y-direction moving distance |
|--------------|-----------------------------|-----------------------------|
| 165 x 210 mm | 60 mm                       | 30 mm                       |

Double Lens Condenser: Numeric opening = 1,2 with built-in Iris cover  
Working distance of stage (vertical) 18 mm  
Range of fine focusing: 14 mm, the graduation is 0,002 mm.  
Range of coarse focusing: 14 mm

### b) Components of Student Microscope M-TOP 600 and M-TOP 600 LED

Study the picture below to become familiar with the different parts of your microscope.

### c) Equipment

- 1 Dust Cover
- 1 Color Filter
- 1 Spare Halogen Bulb 20 Watt, 6 Volt for M-TOP 600
- 1 Battery Recharger for M-TOP 600 LED

### Intended use of the microscope

M-TOP 600 microscopes are produced for clinical routine check-ups as well as for teaching purposes. They are the ideal choice for biological, bacteriological, pathological, and pharmaceutical research. These microscopes are equipped with achromatic objectives, eyepieces, Abbe condensers, rough and fine adjustment. The model M-TOP 600 LED additionally features a new, improved 4,5 V LED illumination system, which enables a more softer and convenient top view.

## Directions for use

### a) Preparations: Setting up the microscope

Take the microscope out of its Styrofoam wrapping and put it onto a stable and level working surface.

### b) Work Instructions

#### Illumination

The model M-TOP 600 works with a halogen bulb (20 Watt, 6 Volt) imbedded in the base. This bulb has an average operating time of 3000 working hours.

The model M-TOP 600 LED is equipped with LED lights. The electricity for the lighting is stored in rechargeable batteries. The light of the LED bulbs hardly heats liquid specimens. This bulb has an average operating time of 100 000 working hours. The lighting is switched on at the back of the base. The rechargeable batteries are best recharged overnight.

#### Power supply

The model M-TOP 600 LED can be connected to a normal power socket or used cable-free with the rechargeable batteries. One battery charge lasts for about 10 operating hours. The included recharger can only be connected to a power circuit of 220 Volt. The recharger is protected against overcharging. First connect the recharger to the microscope, and then connect it to the power circuit. Should the rechargeable batteries totally be discharged, please recharge them overnight. The batteries can only be recharged when the microscope is switched off. Please do not exchange defective accumulator batteries by yourself. Contact our service personnel listed under "Guarantee".

The model M-TOP 600 only operates with a power network of 220 Volt.

#### Placing the object plate

Place the object plate provided with an object sample onto the object stage, and fasten it carefully with the help of the mobile spring clips. By turning the y/y object stage adjustment knobs (at the left under the object stage), you can move the object to be viewed into the centre of the object stage. Please note that these knobs have two functions: They move the object plate on the object stage forwards and backwards, or from the right to the left.

#### Adjusting the image sharpness

Start by viewing your object with a four times magnification. Turn the button for coarse adjusting until you receive a clear picture. Then use the button for fine adjusting to achieve an optimal sharpness. Proceed the same way with every other object. Pay attention that the object does not touch the objective lens.

#### Brightness control

There are several ways of adjusting the illumination projected onto the object. Light influences the viewing quality especially with high magnification. To regulate the illumination LED through the condenser, the blind can be opened or shut. To move the light beam closer to or farther from the object, the condenser part can be turned up or down with the regulating knob. You can also adjust the illumination with the dimmer. Adjust the dimmer to its lowest intensity, switch on the current, and increase the illumination up to the wanted intensity.

Use these different possibilities of light adjustment in order to reach a perfect top view of your object.

#### Changing the objective

When you reached the optimal top view with the smallest magnification (4x), turn the objective lens piece to the next magnification (10x). The objective lens piece should click into position. Only a slight adjustment should be necessary, to receive a clear picture again. Should you increase the magnification further; the top view will be improved by adjusting the light.

### **Use of the object stage**

The Betzold microscopes are provided with a lock for the object stage. This allows only a certain maximum height of the stage to avoid collisions with the objectives. This lock is already adjusted correctly, and should not need any further adjustment. However, should you have the impression that you need to approach the objective nearer, you can adjust the lock yourself. Please pay attention when doing this as damaged objectives will not be replaced, as the damage has been CAUSED BY THE CUSTOMER! So, always keep in mind to not screw up the object plate too high!  
When using 60x magnification, the lens appears to be very close the slide holder.

### **After use**

When having finished working with the microscope, turn off the current. Use the button for coarse adjusting to move the object stage downwards. The device should be protected by a dust cover, or disassembled and put back into the original Styrofoam packing.

### **c) Exchanging the halogen bulb of the microscope Model M-TOP 600**

A halogen bulb 20 Watt, 6 Volt is required.

First pull out the power plug! Loosen the screws at the bottom of the microscope, and remove the base plate. Replace the used bulb by a new one. Pay attention – the bulb should not be touched with your fingers, as this might shorten the operating time. Replace the base plate.

### **Suggestions for the observation of specimen**

Children should not experiment without adult supervision!

1. Place some grains of salt or sugar on the object holder.
2. Dissolve a teaspoon of salt/sugar in warm water in a test tube. Then use the pipette to place some drops on an object holder. Let them dry and observe the recrystallization of the salt/sugar.
3. Shake some pollen of a blossom onto an empty object holder.
4. It is interesting to observe parts of feathers, or any hair, that of domestic animals included. Try to compare the hair of different animals. You can fix the hair with a small piece of sticky tape onto the object holder.
5. Observe the fibres of different types of paper, or compare the print of a newspaper with that of a type-writer or a laser-printer.
6. Look at the colours in comics, or at strings and fibres of different materials.
7. Compare natural fibres like cotton wool or wool, to synthetic ones like polyester.
8. Water, out of a pond, contains living organisms. Put a drop of it onto the object holder with the dip to observe them.
9. Very small insects can be looked at with little magnification.
10. Or try to look at parts of bigger beetles, like legs or wings. Doing this, the needle is useful for the positioning. Always take insects with tweezers, and wash your hands carefully after handling them.
11. Watch the mould on old bread or other food. Or create your own mould culture, by putting a little piece of bread into a petri dish for some days. A little moistness accelerates the moulding process. (Here too, wash your hands carefully after work.)
12. Many objects seem to be transparent under the microscope. Usually these objects are dyed, to make the cells more visible. A current dye is iodine. You can also dip your tissue sample for some minutes into normal, diluted food colouring, in order to dye it. Use a pair of tweezers to transport the dyed objects. Bear in mind that the colour will dye any objects contacted. Therefore take care of clothing and furniture.

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